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healthcare professionals, and conducted in laboratories, GP offices, clinics, pharmacies, etc.

7. The Test Cartridge should remain in its original sealed pouch until ready to use. Do not use the Test Cartridge if the pouch is punctured or not well sealed. Discard after single use.
8. There is a blue line on the test membrane, it will disappear after sample adding. This indicates that Test Cartridge has been used. Do not reuse the Test Cartridge.
9. Do not use damaged or stained materials provided in the test kit.
10. The test kit and instrument should be used away from vibration and magnetic field. During normal usage, the instrument may introduce minute vibration, which should be regarded normal.
11. The Pipette Tips and Detection Buffer Tubes should be used for one specimen only. Discard after single use.
12. Do not use whole blood specimen when hemolysis or blood clot appears.
13. Do not smoke, eat, or drink in areas in which specimens or kit reagents are handled.
14. The blood specimens and used materials, such as Test Cartridges, Detection Buffer Tubes and Pipette Tips, are potentially infectious. Proper safety techniques, handling and disposal methods should be followed in accordance with standard procedures and relevant regulations observed by biological hazard materials.
15. The Fineware™ Cortisol Rapid Quantitative Test should not be used as absolute evidence for adrenal function. The results should be interpreted by the physician along with clinical findings and other laboratory test results.
16. The test will be applied on a routine basis and not in emergency situations.
17. If you have any questions or need help, please contact the local distributor to solve problems timely.
18. Notice to the users: Any serious incident that has occurred in relation to the Fineware™ Cortisol Rapid Quantitative Test shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

MATERIAL**Material Provided**

1. 25 individual sealed pouches, each containing:
 - Test Cartridge
 - Desiccant Pouch
2. ID Chip
3. Leaflet with Instructions For Use
4. 25 Pipette Tips (for 100 µL transfer pipette set)
5. 25 Detection Buffer Tubes

Material Required But Not Provided

Function	Material Name	Note
Test instrument	Fineware™ FIA Meter	Model No.: FS-112
	Fineware™ FIA Meter Plus	Model No.: FS-113
	Fineware™ FIA Meter II Plus SE	Model No.: FS-114
	Fineware™ FIA Meter III	Model No.: FS-205

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about their accuracy.

Note: Please refer to the Instructions For Use of Finecare™ Cortisol Control for detailed operation.

TRACEABILITY

Finecare™ Cortisol Rapid Quantitative Test has been standardized against the IRMM/IFCC-451 Panel (ID/GC/MS)

LIMITATIONS OF PROCEDURE

1. This test has been developed for testing human whole blood, serum and plasma specimen only.
2. This test has been verified for healthcare professionals use. Refer to point 6 of the PRECAUTIONS in this instruction for requirements of training and qualifications required by the users. Note that this test is not used for self-testing.
3. The results of Finecare™ Cortisol Rapid Quantitative Test should be evaluated with all clinical and laboratory data available. If the test results do not agree with the clinical evaluation, additional tests should be performed.
4. The false positive results include non-specific adhesion of some components in specimen that have similar epitopes to bind captured and detector antibodies.
5. In the case of false negative results, the most common factors are: non-responsiveness of antigen to the antibodies by that certain unknown components are masking its epitope, such that antigen cannot be seen by the antibodies; instability of Cortisol, resulting in degradation with time or temperature, such that they become no longer recognizable by antibodies.
6. Whole blood using anticoagulants other than EDTA, Heparin Sodium or Sodium Citrate has not been evaluated in Finecare™ Cortisol Rapid Quantitative Test and thus should not be used.
7. Other factors may interfere with Finecare™ Cortisol Rapid Quantitative Test and may cause erroneous results. These include technical or procedural errors, as well as additional substances in blood specimens.

PERFORMANCE CHARACTERISTICS

Accuracy

A comparative study was tested for 108 clinical samples using Finecare™ Cortisol Rapid Quantitative Test and the Roche Cortisol assay. The Correlation Coefficient (r) was 0.986.

Measuring Range and Detection Capability

- **Measuring Range:** 50–1000 nmo/L
- **Limit of Detection (LoD):** 50 nmol/L

Precision

Intra-Lot Precision:

Within-lot precision has been determined by Cortisol precision controls with

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INDEX OF SYMBOLS

	See Instructions for Use		Tests per kit		Manufacturing date
--	--------------------------	--	---------------	--	--------------------

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Guangzhou Wondfo Biotech Co., Ltd.