



深圳市新产业生物医学工程股份有限公司
Shenzhen New Industries Biomedical Engineering Co., Ltd(Snibe)

Declaration

To whom it may concern,

We, Shenzhen New Industries Biomedical Engineering Co., Ltd., certify that the MAGLUMI Free PSA(CLIA) (REF:130201035M, 130601035M, 130701035M) can use the reference range as follow.

Males: < 1.5 ng/mL(95th percentile)

Females: < 0.1 ng/mL(95th percentile)

Shenzhen New Industries Biomedical Engineering Co., Ltd.



2024-06-13